


2025

Product Catalogue



Contents

pg. 2	About us & Contacts
pg. 3	Handheld Ultrasound Devices
pg. 4	Medical Paper ◦ Recording Paper ◦ Thermal Ultrasound paper
pg. 5	ECG Electrodes
pg. 7	Defib & Neutral Electrodes ◦ Neutral Electrodes ◦ Defibrillation Electrodes
pg. 8	Sensors & Accessories ◦ BP Cuffs ◦ SP02 Sensors ◦ Pressure Infusion Cuffs
pg. 9	Airway ◦ Oral Care ◦ Nasal & Laryngo-Tracheal Atomisers
pg. 10-11	Tracheostomy Kits ◦ Ciaglia/Dilator Technique Kits ◦ Griggs Technique Kits ◦ Tube Exchange Guide ◦ Basic Tracheostomy Kits ◦ Kit Contents
pg. 11-12	Electrosurgery ◦ Smoke Evacuation ◦ Diathermy Pencils ◦ Blade, Needle, Loop, & Ball Electrodes
pg. 13	Our Brands & Suppliers


About Us

Anchored in excellence, Driven by care

With over a decade of trusted service, we are ever committed to building trust and supporting our customers. Customers are our core focus, and our responsive team is dedicated to understanding and meeting your unique needs. Backed by a robust distribution network, Anchor Medical ensures efficient service to all healthcare facilities across South Africa.

Quality Assurance

Quality is our top priority when it comes to our products. We meticulously curate our product range, ensuring that we supply only the highest-quality medical devices to our clients. Every product we offer provides absolute reliability and safety. We pride ourselves in delivering solutions that healthcare professionals can trust, offering the optimal care for their patients.



Contact us

Contact Details:

Email: Sales@anchormed.co.za
Website: www.anchormed.co.za
Telephone:
Gauteng Head Office:
(+27) 11 262 5045

Address:

Anchor Medical (PTY) LTD
17 Carey Street,
Wynberg EXT 1,
Sandton, 2090
Johannesburg
South Africa

Handheld Ultrasound Devices

Meet the **Vinno Q+**, Vinno's 3rd generation of handheld ultrasound. 30% smaller and 50% lighter than before, makes it the ideal companion for sideline care. With a user oriented Android App, the Q+ allows for effortless real-time diagnosis wherever you are.

Key features include:

- Windows & Android compatibility
- B Mode, Colour Doppler, M Mode & PW Mode
- Save images to the cloud with the FlyinSono App
- One handed operation with buttons to change gain, frequency and depth
- IPX7 waterproof body for easy cleaning & disinfection
- Weighs only 135 g



Code	Probe Type	Frequency	Max Depth
Vinno Q+ 7L	Linear Probe	6-16 MHz	6 cm
Vinno Q+ 2p	Phased Array	2-5 MHz	36 cm
Vinno Q+ 3C	Convex Probe	2-6.5 MHz	36 cm

Portability without Compromise

The Q+ maintains excellent image quality despite its small size and weight, making it the perfect device for Point-of-Care Ultrasound. Why accept any compromises?



Liver



Thyroid



Four Chamber view



Liver CF



Kidney CF



Carotid CF

Recording & Thermal paper

EasyTrace recording paper offers excellent performance at a fraction of the price. For this reason, we have become the market leader in South Africa for recording paper. Healthcare providers have come to rely on the clear and reliable images our paper provides. We supply ECG, CTG and Defib paper for all common printers.



Recording Paper

Code	Printer	Size
FQW210-3-140	Nihon Kohden Cardiofax	210x140 mm, 215 sheets
M1911A	Philips M1911A	152x90mm, 150 sheets
FQS50-3-100	Nihon Kohden FQS50.3.100	50x100 mm, 300 sheets
4305BA0/DAO	Corometrics 4305BAD/DAO	152x90, 150 sheets
DFM100	Philips DFM100	50mm roll, 20m
M1709A	Philips M1709A/M2483A	210x300 mm, 200 sheets
Zoll 8000	Zoll 8000-03000	80mm roll, 20m
Edan F3/F9	Edan F2/F3/F9	152x90mm, 150 sheets
EDAN SE1200	Edan SE12, MAC1200	210x297mm , 150 sheets
WA 94018	Welch-Allyn CP100	210x280mmx, 200 sheets
40457C/D	HP 40457A/B/C	50mm roll, 30m

more options available
upon request

Thermal papers



Code	Brand	Type
UPP110S/HD/HG	Sony	Thermal Roll, 110mm, various
UPP210HD	Sony	Thermal Roll, 210mm, High Density
KH65M	Mitsubishi	Thermal Roll, 110mm, High Density
UPC-21L	Sony	Print Media pack, A6
UPC-R80MD	Sony	Self Laminating print media pack, A4

ECG Electrodes

Skintact ECG electrodes have been a mainstay in South Africa for the last 20 years. Skintact® ECG's are manufactured in the Austrian Alps using custom built machinery for only the best quality. This way they ensure absolute reliability, ease of use, and comfort for the patient. Our extensive range has electrodes for every application including low-cost multipurpose electrodes, stress test electrodes, and monitoring electrodes.



	Code	Description	Size	Gel	Adhesion	Application
	FS-RG1/10	Multipurpose electrode for adult and paed	41x32 mm	Solid Gel	Strong	Resting, monitoring
	FS-RB4/5	Foam Electrode with Banana plug connector	44x28 mm	Liquid gel	Strong	Diagnostic
	F-521	Foam Electrode, Water repellant	50x35 mm	Solid gel	Strong	Monitoring, ideal for operating room
	F-50C	Foam electrode, Carbon Stud	Ø 50 mm	Liquid gel	Medium	X-Ray, MRI Conditional
	FS-50	Foam electrode, Water repellant	Ø 50 mm	Liquid gel	Medium	Echochardiogram, Emergency room
	F-401	Foam electrode for paediatrics	Ø 40 mm	Solid gel	Medium	Resting, monitoring
	F-301	Foam electrode for neonates	Ø 30 mm	Solid gel	Medium	Resting, Monitoring
	PD50-SFC	Prewired foam electrode for neonatal	Ø 30 mm	Solid gel	Medium	Resting, Monitoring, X-ray

Nuetral & Defibrillation Electrodes



Our **Skintact®** Neutral Electrodes and Defib pads offer complete reliability and safety in high risk situations. Both pads and electrodes contain overlapping gel technology to minimize burns from the exposed aluminium edge. The geometry and gel used ensure excellent adhesion and contact to the patient. Our Neutral electrodes further are designed with round edges to disperse currents more evenly, leading to a cooler interface.

Neutral Electrodes

Code	Type	Area	Weight Class
WR21A30	Split Plate, Corded, Adult & Paed	90 cm ²	>5kg
WR21	Split Plate, Uncorded, Adult & Paed	90 cm ²	>5kg
WR01	Solid Plate, Uncorded, Adult & Paed	90 cm ²	>5kg
RSW23A30	Split Plate, Corded, Neonates	50 cm ²	<5kg



Defibrillation Electrodes



Code	Compatible Device	Size
DF27N	Philips Heartstart: FR2, FR2+, FR3 (no electrode self-test) With adapter cable M3508A or With adapter cable M3507A and adapter 05-10200: Philips Heartstart MRx, XL, XL+ / Efficia DFM100	Adult/Paed
DF28N	Zoll: E Series, M Series, R Series, X Series	Adult/Paed
DF43N	Nihon Kohden: With pad adaptor JC-755V: TEC-5500 Series, TEC-7600 Series, TEC-7700 Series With pad adaptor JC-855V: TEC-8300 Series	Adult/Paed
DF32G	Philips Heartstart: With adapter cable M3508A or With adapter cable M3507A and adapter 05-10200: Philips Heartstart MRx, XL, XL+ / Efficia DFM100	Paed/Neonates
DF37G	Nihon Kohden: With pad adaptor JC-755V: TEC-5500 Series, TEC-7600 Series, TEC-7700 Series With pad adaptor JC-855V: TEC-8300 Series	Paed/Neonates

Sensors & Accessories

FMT offers a wide range of SP02 sensors, BP Cuffs, & Pressure Infusion cuffs. All products are manufactured to the highest quality ensuring adherence to the highest of standards. Each range includes reusable & disposable/single-patient options with multiple adapters & connectors, ensuring they work in every application.



BP Cuffs



Code	Type
FMT-S/FMT-D	Reusable Cuff with Bladder, Single/Double Lumen
FMT-SB/FMT-DB	Reusable Cuff w/o Bladder, Single/Double Lumen
FMT-DS/FMT-DD*	Single Patient Cuff w/o Bladder, Single/Double Lumen

Multiple connectors & sizes available for each.

*Also available as Neonatal BP Cuffs.

SpO2 Sensor

Code	Type
FMT-RAF/B	Adult Finger Clip Sensor
FMT-RPF	Paediatric Finger Clip sensor
FMT-RAS/B	Adult Soft Finger Sensor
FMT-RPS	Paediatric Soft Finger Sensor
FMT-RIS	Infant Soft Finger Sensor
FMT-RNS	Neonatal Soft Finger Sensor



Multiple Connectors & Sizes Available for each.

Pressure Infusion Cuffs



Code	Type
FMT-PI-05	Reusable Pressure Infusor, 500 ml
FMT-PI-10	Reusable Pressure Infusor, 1000 ml
FMT-PI-30	Reusable Pressure Infusor, 3000 ml
FMT- DPI / 05	Disposable Pressure Infusor, 500 ml
FMT- DPI / 10	Disposable Pressure Infusor, 1000 ml

Airway

The **Atorex** atomisation range offers a safe, easy, and fast method for delivering medications and topical Anaesthesia via the airway. It reduces the need for intravenous access that can often be painful and upsetting for the patient. Delivery of atomised medication is also relatively painless and requires minimal caregiver training.



Nasal & Laryngo Tracheal Atomisers



Code	Type
AN01	Nasal Atomiser
AN02 (V)	Nasal Atomiser with 3ml syringe (with optional vial adapter)
AN03 (V)	Nasal Atomiser with 1ml syringe (with optional vial adapter)
AT01-11/21	Laryngo Tracheal Atomiser, 11cm/21cm extension
AT01-11/21	Laryngo Tracheal Atomiser, 11cm/21cm extension
AT02-11/21 (V)	Laryngo Tracheal Atomiser with 3ml syringe, 11cm/21cm extension (with optional vial adapter)

Oral Care Devices

The Deas Oral Care Line offers a complete solution for the prevention of Ventilator Acquired Pneumonia (VAP), with the added benefit of improved oral health & patient comfort.

Code	Type
SS01	Oral Care Swab
SA01 NS	Oral Care Swab with Suction
09488 NS	Oral Care kit with Oral Swab & Suction Toothbrush



PDT Kits

The **Accex®** range offers a comprehensive range of kits for percutaneous tracheostomy at the bedside or in the OR. The range includes both Ciaglia (Dilator) and Griggs technique kits:



ACCEx
TRACHEOSTOMY KITS

- The Ciaglia kits offer a simple and effortless solution with its grooved dilator to enhance operator's control
- The Griggs kits cost-effective and rapid solutions with its unique dilating forceps allowing for a single step dilation.

Ciaglia Technique



Code	Type	Size
TP01-00	Kit with tracheostomy tube introducers	N/A
TP01-01	Kit w/o tracheostomy tube introducers	N/A
TP90-07/08/09	Kit with standard tracheostomy tube	7 mm, 8 mm, 9 mm
TP90-17/18/19	Kit with reinforced tracheostomy tube	7 mm, 8 mm, 9 mm
TP90-27/28/29	Kit with fixed flanged tracheostomy tube	7 mm, 8 mm, 9 mm
TP90-37/38/39	Kit with fixed flange tracheostomy tube for subglottic suctioning	7 mm, 8 mm, 9 mm

Griggs Technique

Code	Type	Size
TP02-27/28/29	Kit with standard tracheostomy tube	7 mm, 8 mm, 9 mm
TP02-37/38/39	Kit with reinforced tracheostomy tube	7 mm, 8 mm, 9 mm
TP02-07/08/09	Kit with standard tracheostomy tube & dilating forceps	7 mm, 8 mm, 9 mm
TP02-17/18/19	Kit with reinforced tracheostomy tube & dilating forceps	7 mm, 8 mm, 9 mm
TP02-77/78/79	Kit with fixed flanged tracheostomy tube	7 mm, 8 mm, 9 mm
TP02-67/68/69	Kit with fixed flange tracheostomy tube & dilating forceps	7 mm, 8 mm, 9 mm
TP02-74/75/76	Kit with fixed flange tracheostomy tube for subglottic suctioning	7 mm, 8 mm, 9 mm
TP02-64/65/66	Kit with fixed flange tracheostomy tube for subglottic suctioning & dilating forceps	7 mm, 8 mm, 9 mm



PDT Kits

The Accex Tracheostomy Tube Exchange Guide is designed for the safe replacement of tracheostomy tubes, using the seldinger technique with the ability to ventilate/oxygenate during the procedure



Removable
15M connector for
ventilation/oxygenation



Atraumatic rounded closed
tip with holes



Depth markings
with centimeter scale

Tube Exchange Guide

Code	Dimensions	Length
GSC06400	4.0 mm I.D, 6.0 mm O.D	400 mm
GSC07400	4.7 mm I.D, 7.0 mm O.D	400 mm

Basic Tracheostomy Kits

Code	Type	Size
TP20-07/08/09	Kit with standard adjustable flange tracheostomy tube	7 mm, 8 mm, 9 mm
TP21-07/08/09	Kit with reinforced adjustable flange tracheostomy tube	7 mm, 8 mm, 9 mm
TP20-07/08/09	Kit with standard fixed flange tracheostomy tube	7 mm, 8 mm, 9 mm
TP21-07/08/09	Kit with fixed flange tracheostomy tube for subglottic suctioning	7 mm, 8 mm, 9 mm



Kit Contents

Basic Kit Contents:

- Tracheostomy Tube (Variable)
- Tracheostomy Tuber Introducer
- Guiding Catheter
- Adjustable fixation band with Velcro® fasteners
- 9 x 9 cm tracheal compress
- Sterile lubricating jelly sachet

Advanced Kit Contents:

- All contents of basic Kit
- Drape with clear plastic window
- Instrument table cover
- Scalpel no. 15
- 15 G introducer needle
- 16 G sheathed needle
- 14 Fr dilator
- 0.05" J-tip guidewire
- 10 x 10 cm gauzes
- 9 x 9 cm tracheal compress

Electrosurgery

Cimpax offers a wide range of the highest quality electrosurgery electrodes & smoke evacuation handpieces. The C-Vac handpieces provide exceptional performance and comfort for the user, resulting in a safer procedures not only for the patient, but the theatre staff as well. All electrodes are coated in advanced teflon or silicon coatings, significantly reducing time spent cleaning the electrode.

C I M P A X 



30

One day's exposure to surgical smoke is equivalent to smoking 27 to 30 unfiltered cigarettes



150


Hazardous chemicals are found in surgical smoke, including benzene and formaldehyde



2x


Perioperative nurses have twice the incidence of many respiratory problems as compared to the general population

Smoke Evacuation



Code	Description
C31-19-101	C-Vac Tornado Smoke evacuation pencil, coated
C31-19-201	C-Vac Tornado Smoke evacuation pencil, Silicone
C31-20-101	Telescopic smoke evacuation pencil, coated

Diathermy Pencils



Code	Description
PN03-080	Diathermy pencil, coated, 3m cable
PN03-080	Diathermy pencil, uncoated, 3m cable
TC1001	Scratch pad, 5x5 cm

Electrosurgery

Blade Electrodes



C31-10-061 Coated blade electrode 65mm



C31-10-063 Silicone blade electrode 65mm



C31-10-101 Coated blade electrode 100mm



C31-10-161 Coated blade electrode 165mm

(more options available upon request)

Ball Electrodes



C31-12-061 Coated ball electrode 5mm/55mm



C31-12-121 Coated ball electrode 5mm/120mm

Needle Electrodes



C31-11-061 Coated needle electrode 65mm



C31-11-063 Micro-dissection (tungsten) needle electrode 65mm



C31-11-101 Coated needle electrode 100mm

(more options available upon request)

Loop Electrodes



C31-17-120
Loop electrode 10x10mm/120mm



C31-17-121
Loop electrode 15x12mm/120mm



C31-17-122
Loop electrode 20x15mm/120mm

(All loop electrodes colour coded with holder)

Our Suppliers & Brands



DEAS



ACCEX



EasyTrace®

C I M P A X 

ATOREX®



FMT
Medical Sensors, Cables and Accessories



SONY