

Nuetral & Defibrillation Electrodes



Our **Skintact®** Neutral Electrodes and Defib pads offer complete reliability and safety in high risk situations. Both pads and electrodes contain overlapping gel technology to minimize burns from the exposed aluminium edge. The geometry and gel used ensure excellent adhesion and contact to the patient. Our Neutral electrodes further are designed with round edges to disperse currents more evenly, leading to a cooler interface.

Neutral Electrodes

Code	Type	Area	Weight Class
WR21A30	Split Plate, Corded, Adult & Paed	90 cm ²	>5kg
WR21	Split Plate, Uncorded, Adult & Paed	90 cm ²	>5kg
WR01	Solid Plate, Uncorded, Adult & Paed	90 cm ²	>5kg
RSW23A30	Split Plate, Corded, Neonates	50 cm ²	<5kg



Defibrillation Electrodes



Code	Compatible Device	Size
DF27N	Philips Heartstart: FR2, FR2+, FR3 (no electrode self-test) With adapter cable M3508A or With adapter cable M3507A and adapter 05-10200: Philips Heartstart MRx, XL, XL+ / Efficia DFM100	Adult/Paed
DF28N	Zoll: E Series, M Series, R Series, X Series	Adult/Paed
DF43N	Nihon Kohden: With pad adaptor JC-755V: TEC-5500 Series, TEC-7600 Series, TEC-7700 Series With pad adaptor JC-855V: TEC-8300 Series	Adult/Paed
DF32G	Philips Heartstart: With adapter cable M3508A or With adapter cable M3507A and adapter 05-10200: Philips Heartstart MRx, XL, XL+ / Efficia DFM100	Paed/Neonates
DF37G	Nihon Kohden: With pad adaptor JC-755V: TEC-5500 Series, TEC-7600 Series, TEC-7700 Series With pad adaptor JC-855V: TEC-8300 Series	Paed/Neonates