



ECG Electrodes

Skintact ECG electrodes have been a mainstay in South Africa for the last 20 years. Skintact® ECG's are manufactured in the Austrian Alps using custom built machinery for only the best quality. This way they ensure absolute reliability, ease of use, and comfort for the patient. Our extensive range has electrodes for every application including low-cost multipurpose electrodes, stress test electrodes, and monitoring electrodes.



	Code	Description	Size	Gel	Adhesion	Application
	FS-RG1/10	Multipurpose electrode for adult and paed	41x32 mm	Solid Gel	Strong	Resting, monitoring
	FS-RB4/5	Foam Electrode with Banana plug connector	44x28 mm	Liquid gel	Strong	Diagnostic
	F-521	Foam Electrode, Water repellant	50x35 mm	Solid gel	Strong	Monitoring, ideal for operating room
	F-50C	Foam electrode, Carbon Stud	Ø 50 mm	Liquid gel	Medium	X-Ray, MRI Conditional
	FS-50	Foam electrode, Water repellant	Ø 50 mm	Liquid gel	Medium	Echochardiogram, Emergency room
	F-401	Foam electrode for paediatrics	Ø 40 mm	Solid gel	Medium	Resting, monitoring
	F-301	Foam electrode for neonates	Ø 30 mm	Solid gel	Medium	Resting, Monitoring
	PD50-SFC	Prewired foam electrode for neonatal	Ø 30 mm	Solid gel	Medium	Resting, Monitoring, X-ray